

Serological Pipettes

NEST serological pipettes can be used to liquid transfer purpose of both tissue culture and bacterial culture experiments.

Features

1. Made of high clarity polystyrene.
2. Sterilized by E-beam.
3. Non-pyrogenic, non-cytotoxic and non-hemolytic.
4. DNase/RNase free.
5. Polyolefin fiber filter to reduce contamination, meanwhile it can help prevent overflow.
6. Edged color-coded band for identification.
7. Exceptionally crisp, dark, easy to read printed graduation.
8. Reverse graduations, negative graduations for extra capacity.

Parameter

Volume	Graduations	Negative Graduation	Color Code
1 mL	1/100 mL	0.4 mL	Yellow
2 mL	1/50 mL	0.6 mL	Green
5 mL	1/10 mL	2.5 mL	Blue
10 mL	1/10 mL	3 mL	Orange
25 mL	2/10 mL	7 mL	Red
50 mL	5/10 mL	6 mL	Black
100 mL	1 mL	0 mL	Pink



Product Range



Plastic/Paper Individually Wrapped, Sterile

Cat. No.	Description	/Pk	/Cs
324001	1mL Serological Pipette	500	3000
325001	2mL Serological Pipette	400	2400
325011	2mL Aspirating Pipette, No Printing, No plug	200	800
326001	5mL Serological Pipette	200	800
327001	10mL Serological Pipette	200	800
328001	25mL Serological Pipette	200	800
329001	50mL Serological Pipette	100	600
329501	100mL Serological Pipette	50	300




Bulk, 20 per bulk bag, Sterile

Cat. No.	Description	/Pk	/Cs
324002	1mL Serological Pipette	1000	3000
325002	2mL Serological Pipette	1000	3000
326002	5mL Serological Pipette	500	1500
327002	10mL Serological Pipette	500	1500
328002	25mL Serological Pipette	200	800
329002	50mL Serological Pipette	100	600

Technical Drawing



NEST China	Tel : 0510-6878 8698 / 6878 8718	NEST USA	Tel : +1-732-381-0268
	Email : info@nest-wuxi.com		Email : sales@nestscientificusa.com
	Website : www.cell-nest.com		Website : www.nestscientificusa.com
	Add : No.530, xida road, new district, Wuxi, Jiangsu, China		Add : 3 Convery Blvd, Woodbridge, NJ 07095 USA



NEST China	Tel : 0510-6878 8698 / 6878 8718	NEST USA	Tel : +1-732-381-0268
	Email : info@nest-wuxi.com		Email : sales@nestscientificusa.com
	Website : www.cell-nest.com		Website : www.nestscientificusa.com
	Add : No.530, xida road, new district, Wuxi, Jiangsu, China		Add : 3 Convery Blvd, Woodbridge, NJ 07095 USA