

# LTS Pipette Tips

## Introduction

With the constant efforts and many tests of our engineers, we are glad to announce that the alternative RAININ LTS tips are also available for purchase.

This kind of tip is different from the universal tips. With the similar design to RAININ LTS Pipette Tips, it could be compatible well with RAININ LTS pipettes.

Meanwhile, like all other pipette tips, it is made from virgin polypropylene and manufactured in Class 100,000 clean room. Also, our smart production line ensures minimal human contact during the production. Thus, these tips are free of DNase, RNase, human DNA, and endotoxins.

## Features

1. Made from virgin polypropylene.
2. Compatible well with RAININ LTS pipettes.
3. RNase, DNase, Human DNA, Pyrogen-Free.
4. Sterilized by E-beam SAL=10<sup>-6</sup>

## Product Description

### Rainin LTS Tips - Filtered, Sterile

Cat.No.	Description	/pk	/cs
311061	Rainin LTS 10uL Filter Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800
312061	Rainin LTS 200uL Filter Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800
319061	Rainin LTS 300uL Filter Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800
313061	Rainin LTS 1000uL Filter Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800

**NEST  
China**

Tel : 0510-6878 8698 / 6878 8718  
 Email : info@nest-wuxi.com  
 Website : www.cell-nest.com  
 Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST  
USA**


Tel : +1-732-381-0268  
 Email : sales@nestscientificusa.com  
 Website : www.nestscientificusa.com  
 Add : 3 Convery Blvd, Woodbridge, NJ 07095 USA

**Rainin LTS Tips - Non-filtered, Sterile**


Cat.No.	Description	/pk	/cs
301061	Rainin LTS 10uL Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800
302061	Rainin LTS 200uL Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800
305061	Rainin LTS 300uL Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800
303061	Rainin LTS 1000uL Pipette Tip, Clear, Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	960	4800

**Technical Drawing**


**10uL Filter Pipette Tip**




200uL Filter Pipette Tip



300uL Filter Pipette Tip



1000uL Filter Pipette Tip




**NEST  
China**

Tel : 0510-6878 8698 / 6878 8718  
 Email : info@nest-wuxi.com  
 Website : www.cell-nest.com  
 Add : No.530, xida road, new district, Wuxi, Jiangsu, China


**NEST  
USA**

Tel : +1-732-381-0268  
 Email : sales@nestscientificusa.com  
 Website : www.nestsscientificusa.com  
 Add : 3 Convery Blvd, Woodbridge, NJ 07095 USA


10uL Pipette Tip



200uL Pipette Tip



300uL Pipette Tip




**NEST  
China**

Tel : 0510-6878 8698 / 6878 8718  
 Email : info@nest-wuxi.com  
 Website : www.cell-nest.com  
 Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST  
USA**

Tel : +1-732-381-0268  
 Email : sales@nestscientificusa.com  
 Website : www.nestscientificusa.com  
 Add : 3 Convery Blvd, Woodbridge, NJ 07095 USA

1000uL Pipette Tip



**NEST  
China**

Tel : 0510-6878 8698 / 6878 8718  
 Email : info@nest-wuxi.com  
 Website : www.cell-nest.com  
 Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST  
USA**

Tel : +1-732-381-0268  
 Email : sales@nestscientificusa.com  
 Website : www.nestscientificusa.com  
 Add : 3 Convery Blvd, Woodbridge, NJ 07095 USA